




Compatible systems

- | | |
|-------------------------|----------------------------|
| // 3i Certain | // Straumann bone level |
| // Frialit | // Astratech |
| // Nobel biocare active | // 3i Ossetite |
| // Straumann ITI | // Nobel biocare replace |
| // Zimmer | // Nobel biocare branemark |



Dental House, 239 Staines Road West, Sunbury-on-Thames,
Middlesex TW16 7BH Tel. 01932 786 296
chris@chrishowe.co.uk www.chrishoweceramics.co.uk



chameleon[®]
implants
strength in your
surroundings



Are you looking for the best aesthetic solution for cement and screw retained implant prosthetics? Chameleon Zirconia custom abutments could be the answer.

By utilising our unique virtual abutment design software, the zirconia abutments are individually designed within the final tooth shape. The result is an abutment with outstanding function and aesthetics.

The restorative procedure is simple with reduced chair time, it will also eliminate the need for inventory of stock abutments and have a fixed cost custom abutment since no precious metals are used.

The CAD construction of Chameleon abutments allows a greater accuracy when designing emergence angles and paths of insertion and with a highly polished zirconia collar, provides better acceptance to the soft tissues the result the optimum aesthetic abutment for single and multiple cases.

Chameleon abutments are platform independent and available in 8 vita shades to complement the restoration.

The Chameleon abutment in zirconia is perfect for challenging cases in the aesthetic zone and is strong enough to be used in all positions in the mouth.



Five important factors to consider when choosing Chameleon Zirconia Abutments:

// **Freedom**

Experience the freedom of unlimited possibilities for cement-retained implant prosthetics. Chameleon abutments are available for most implant systems.

// **Aesthetics**

Design, support, function and colour make each Chameleon abutment as individual as your patient. Each abutment is designed to be anatomically correct. The goal is to provide the optimal aesthetic result.

// **Simplicity**

Send fixture level impressions to Terec without spending valuable time on abutment selection or ordering parts. Chameleon abutments eliminate the need for chair side modification and can correct variations in implant angle and depth.

// **Compatibility**

The Chameleon zirconia abutment with its highly polished custom fit collar has an improved acceptance to soft tissue creating a more natural healthy and aesthetic look.

// **Value**

As well as savings in valuable chair side time, Terec also provide a fully customised Chameleon zirconia abutment complete with screw and finished zirconia crown.